

mesotherapy collection

2%hy-ac low viscosity solution

Imunohealth S.L.
Escultor Llimona, 5
E-08328 Alella
+34 930 153 890
www.mesoinstitute.es (English)
sales@mesoinstitute.es
www.mesoinstitute.eu (Spanish)
info@mesoinstitute.eu

Registered in European Union as Cosmetic Products
according to the Cosmetic Products Regulation (EC)
No 1223/2009.
Manufactured according to Good Manufacturing
Practices (GMP).
Made in the European Union.
See the applicable law in your Country for traditional
mesotherapy treatments.

DESCRIPTION

Dehydration and loss of firmness are signs of skin aging that lead to wrinkle formation and sagging of some face areas.

2 % hy-ac low density solution^{/28} is a non cross-linked hyaluronic acid of low molecular weights that acts as anti-aging treatment used to improve the texture of the skin, the fine wrinkles and scars in the face and neck area when is distributed into the mesoderma as well as facilitating the delivery of other substances to the skin.

It's low viscosity allows a better distribution through the dermal layer improving the texture and hydration of the skin avoiding water loss and at the same time stimulating collagen I synthesis.

The skin is therefore progressively tightened and moisturized. The skin is younger, nicer, both flexible and firm.

ACTIVE INGREDIENTS

Hyaluronic acid (combination of HA of high and low molecular weights).

BENEFITS

Reduces wrinkles and fine lines
Stimulates collagen synthesis
Firms the skin
Moisturizes
Improves skin softness

POSSIBLE MIXTURES

vit-F10 cocktail^{/18}
vitamin C^{/32}
detox solution^{/07}
anti-OX solution^{/08}
anti-ROS solution^{/09}

MESO PROTOCOL

Recommended time between sessions

Using meso-pen or derma-roller 7-10 days interval
Traditional 7-10 days interval

Number of sessions

Starting treatment 4-6
Maintenance every 3-4 months



PERESENTATION & REFERENCE

Vial.....5 mL (0.176 fl. oz)
Box.....10 x 5 mL vials
Reference.....102087

Biological product sterilized by filtration
No paraben, alcohol, ingredients of animals origin or colouring.
Microbiologically tested.